

Tired of looking tired?

Eyelid & Brow Lift

Cheek & Jaw Lift

Chin Lift

Neck Lift

Turn back time and wake up to a **new you!**



skinmed^{LTH}
A three dimensional alternative to surgery

The latest non-surgical facelift treatment combining ultrasound (HIFU) and radio frequency to lift and tighten loose or sagging skin on the neck and face ***without any down time***

Ask staff for details

www.3d-lipo.com

@3Dlipo

3D-lipo

The Latest Non Surgical Facelift Treatment

ABOUT

3D-skinmed represents a new Generation of non-surgical, non invasive procedures for the face that utilizes the science of Ultrasound (HIFU) and Radio Frequency and the bodys natural healing process to lift, tone and tighten loose skin. The treatment is natural with no down time and delivers noticeable results

THE SCIENCE

The 3D-skinmed treatment uses the safe time tested energy of Ultrasound (HIFU) which focuses energy to stimulate the deep support layers of the skin without causing damage to the skins surface. The treatment stimulates the growth of new collagen resulting in lifted, firmer and toned appearance over time

WHAT TO EXPECT

Results will develop over the course of 2-3 months with some clients experiencing continued improvement for up to 6 months. Whilst our procedure will not be as drastic as surgical intervention it is an inviting alternative for those who are not ready for surgery.

FREQUENTLY ASKED QUESTIONS.

What is unique about the 3D-skinmed treatment?

3D-skinmed is the only treatment that can specifically target the deep foundation below the skin similar to that of cosmetic surgery without cutting or disrupting the surface of the skin utilizing focused Ultrasound (HIFU). In addition we have adopted a 3 Dimensional approach by utilizing Radio frequency for Superficial skin tightening and Impact technology for deep tissue hydration.

Is 3D-skinmed surgery?

No. Its focussed on the skin and the support layer below it without cutting or damaging the surface of the skin

Can 3D-skinmed replace a face lift?

Although 3D-skinmed is Uplifting it is not a facelift. . Whilst our procedure will not be as drastic as surgical intervention it is an inviting alternative for those who are not ready for surgery.

Is 3D-skinmed safe?

Thousands of treatments have been performed worldwide using this type of technology without any significant adverse events. In addition, ultrasound energy has a proven track record, with use in the field of medicine for more than 50 years

Who is a good candidate for 3D-skinmed?

A good candidate for 3D-skinmed is someone with skin that has relaxed to the point of looking, and often feeling, less firm i.e. tired. Typically, those in their thirties and older who have mild to moderate skin laxity are good candidates. If you are unsure about your compatibility please ask staff for a free consultation.

How long does the treatment take?

A full face and neck treatment takes approximately 60 minutes, whilst a partial face treatment can take around 30 minutes. There is no downtime post treatment.

What problem areas does 3D-skinmed target?

3D-skinmed addresses any area where there are mild to moderate skin laxity issues. Specifically, by reconditioning and renewing the skin from within, 3D-skintech can lift the brow, which in turn reduces excess skin on the lids, opens up the eyes, and gives a more rejuvenated look overall. It is also highly effective on loose skin on the neck area and brings definition to the face and jaw line.

How long until I see results of the 3D-skinmed treatment?

Clients may witness some initial effect but the ultimate lifting and firming will take place over 2-3 months, as collagen is renewed and replaced with new stronger collagen. As this collagen building process continues, continued improvements can appear up to 6 months following a treatment.

How long do the 3D-skinmed results last?

The firming benefits achieved following the treatment should last for over 1 year, but skin continues to age. Yearly touch-up treatments are recommended to combat the ageing process.

How many treatments will I need?

The treatment protocol is on 3D-skinmed Focused Ultrasound (HIFU) treatment followed by 6-8 Radio Frequency treatments at weekly intervals following the initial treatment. For the majority of clients this is sufficient however a further Ultrasound (HIFU) treatment may be required dependant on the extent of the skin laxity and the body's response to treatment

Is there any downtime with 3D-skinmed treatment?

With the 3D-skinmed treatment there is no downtime. You are able to return to your normal activities right away and there are no special aftercare that you need to take.

Will I feel the ultrasound during the treatment?

Yes. What you are feeling during the treatment is actually tiny amounts of concentrated ultrasound energy being delivered below the surface of the skin. Clients have reported that there can be some discomfort while the energy is being applied but it is short lived and a good sign that the collagen restructuring process has been started.

Are there any side effects?

There may be slight redness for about one hour following the procedure, and a low percentage of clients may have slight swelling, tingling or tenderness to touch, but these are mild and temporary in nature.

Ask staff for details

www.3d-lipo.com

@3Dlipo  3D-lipo